

# Liquid Handling



# Single Channel Fixed Volume Pipette

TopPette



## Features

- The low part is available for autoclaving.
- Cover large volumes range from 5µl to 5000µl.
- Long lifetime over 100,000 times pipetting.
- Compatible with most universal brands tips.
- Ergonomic design provides excellent operating experience.
- Online calibration available, easy calibration and maintenance with tool supplied
- High accuracy, error limits are better or equal to international standard ISO 8655-2. Each pipette supplied with individual test certificate.

## Specifications

Volume Range	Test Volume (µl)	Accuracy error		Precision error	
		%	µl	%	µl
5µl	5µl	1.3%	0.065	1.2%	0.06
10µl	10µl	0.8%	0.08	0.8%	0.08
20µl	20µl	0.6%	0.12	0.5%	0.1
25µl	25µl	0.5%	0.125	0.3%	0.075
50µl	50µl	0.5%	0.25	0.3%	0.15
100µl	100µl	0.5%	0.5	0.3%	0.3
200µl	200µl	0.4%	0.8	0.2%	0.4
250µl	250µl	0.4%	1.0	0.2%	0.5
500µl	500µl	0.3%	1.5	0.2%	1.0
1000µl	1000µl	0.3%	3.0	0.2%	2.0
2000µl	2000µl	0.3%	6.0	0.15%	3.0
5000µl	5000µl	0.3%	15	0.15%	7.5