

Liquid Handling



Multi-channel Adjustable Volume Pipette

TopPette

Features



- 8 and 12 channel pipettes are available for standard 96-well plate
- Dispensing head rotates for optimum pipetting convenience
- Individual piston and tip cone assemblies allow easy repair and maintenance
- Compound material tip cone design allows visual seal verification
- Can be used with universal style pipette tips



8 channel

12 channel

Specifications

Picture	Volume Range	Increment	Test Volume (µl)	Accuracy error		Precision error	
				%	µl	%	µl
12-Channel Adjustable Volume Pipette							
	0.5-10µl	0.1µl	10	1.50%	0.15	1.50%	0.15
			5	2.50%	0.125	2.50%	0.125
			1	4.00%	0.04	4.00%	0.04
	5-50µl	0.5µl	50	1.00%	0.5	0.50%	0.25
			25	1.50%	0.0375	1.00%	0.25
			5	3.00%	0.15	2.00%	0.1
	50-300µl	5µl	300	0.70%	2.1	0.25%	0.75
			150	1.00%	1.5	0.50%	0.75
			50	1.50%	0.75	0.80%	0.4
8-Channel Adjustable Volume Pipette							
	0.5-10µl	0.1µl	10	1.50%	0.15	1.50%	0.15
			5	2.50%	0.125	2.50%	0.125
			1	4.00%	0.04	4.00%	0.04
	5-50µl	0.5µl	50	1.00%	0.5	0.50%	0.25
			25	1.50%	0.0375	1.00%	0.25
			5	3.00%	0.03	2.00%	0.1
	50-300µl	5µl	300	0.70%	2.1	0.25%	0.75
			150	1.00%	1.5	0.50%	0.75
			50	1.50%	0.75	0.80%	0.4